

X2

Particle Analyzer CAMSIZER[®] X2

Particle Size and Particle Shape Analysis by
Dynamic Image Analysis



Dry and wet
measurement from
0.8 µm
to **8 mm**


Precise measurement of particle size and shape with the CAMSIZER[®] X2

- Extremely high resolution (0.8 µm/pixel)
- Measuring range from 0.8 µm to 8 mm without hardware adjustment
- Broad dynamic measurement range for broad distributions or multi-modal samples
- Reliable detection of smallest populations (< 0.1 %) of "oversized or undersized" particles
- Particle shape analysis (e. g. for the detection of agglomerates, broken particles, contaminations or rogue particles)
- Innovative optical design provides a resolution of 4.2 megapixel per camera



Particle Analyzer CAMSIZER® X2

Measurement principle	4 – 5
Modular design for optimum measurement conditions	6
Air pressure dispersion Module „X-Dry“ with „X-Jet“	7
Gravity dispersion Module „X-Dry“ with „X-Fall“	8
Wet dispersion Module „X-Flow“	9
Key features	12 – 15
Comparison of measurement techniques	15
Accessories	16
Specifications	16



Experts in Particle Characterization

RETSCH TECHNOLOGY's core competence is to combine innovative particle characterization technology with maximum operating convenience. The instruments, which operate with different measurement techniques, permit analyses of particle size and particle shape for suspensions, emulsions, colloidal systems, powders, granules and bulk materials for quality control and research applications.

The product line covers a size range from 0.3 nm to 30 mm. RETSCH TECHNOLOGY offers tailor-made application and technical services on a global scale.



One measurement method, two measuring ranges

Dynamic Image Analysis with CAMSIZER P4 and CAMSIZER X2

Both CAMSIZER P4 and CAMSIZER X2 are based on the proven and patented dual-camera technology. With more than 1300 units installed worldwide, CAMSIZER systems reliably determine particle size distributions in a range from 1 μm to 30 mm.

While the optics and sample feed system of the CAMSIZER P4 are optimized for the detection of relatively large, free flowing particles (20 μm to 30 mm), the CAMSIZER X2 has been developed for the measurement of fine, agglomerating powders (0.8 μm to 8 mm).

To achieve accurate results for different materials, the CAMSIZER X2 utilizes advanced optical components and offers three particle dispersion options: X-Fall for free fall dispersion (similar to the CAMSIZER P4), X-Jet for air pressure dispersion and X-Flow for liquid dispersion.

Both CAMSIZER systems are equipped with the same powerful and comprehensive software. Which CAMSIZER system is suitable to provide the optimum solution for your analytical requirements, depends on particle size and degree of agglomeration.

For further information visit our website at www.retsch-technology.com/camsizer.



Particle Analyzer

CAMSIZER® X2

Particle Analyzer CAMSIZER X2



Benefits at a glance

- Dynamic image analysis with patented two-camera system (in accordance with ISO 13322-2)
- Wide dynamic measuring range from 0.8 μm to 8 mm
- Newly developed optical system with ultra-strong LEDs for highest resolution and excellent depth of field
- Reliable detection of smallest amounts of "undersized" and "oversized" particles
- Particle shape analysis (e. g. for the detection of agglomerates, broken particles and contaminants)
- Very short measurement time of 1 – 3 minutes
- Excellent reproducibility
- Modular system "X-Change" for dry and wet measurement switching
- Measurement results can be 100% compatible with sieve analysis

The quality control of fine powders can be substantially improved with the CAMSIZER X2: More precise and faster analysis of particle size and particle shape helps to improve product quality, reduce rejected lots and save on costs.

The CAMSIZER X2 is an advancement of the well-proven optical measurement system CAMSIZER for finer samples. The basic differences include higher resolution of the optical systems and improved options for the sample feeding. Fine particles tend to agglomerate which makes it difficult to detect the geometric dimensions of the primary particles. It is therefore beneficial to have various possibilities of feeding the sample into the measurement area. Thus it

is possible to find the best way of dispersing agglomerates without destroying the primary particles for each material. Here, the CAMSIZER X2 provides flexible solutions: from the free fall module, which is the gentlest method for feeding the material, to air pressure dispersion with adjustable pressure and variable nozzle geometry. Additionally, a wet module in which particles are dispersed in various liquids, optionally with ultrasonic dispersion.

Ideal for:

- ▶ **Pharmaceutical powders, granules and fine pellets**
- ▶ **Pulverized and granulated food**
- ▶ **Detergent powder and chemical raw materials**
- ▶ **Plastic powders (even with an electrostatic charge on the surface)**
- ▶ **Metal powders for additive manufacturing, MIM, solder powders**
- ▶ **Abrasives (medium-size and small grit)**
- ▶ **Fine sands and construction materials**
- ▶ **Fine plastic, glass and carbon fibers**

New Optical System


Bright light sources, fast cameras with high resolution, and particularly the new optical split view concept make the CAMSIZER X2 a superior analyzer. The patented split view optical design arranges the optical paths of the two cameras perpendicular to the direction of movement of the dispersed particles. The monochromatic light beams of the two pulsed LED light sources are focused and collimated to illuminate the detection area of each camera respectively (Basic and Zoom). Both light beams are combined with a prism so that they illuminate their respective, overlapping detection areas from the same direction which is orthogonal to the particle beam. This new optical design provides an optimized illumination area, field of view, pulse length and pulse frequency for each camera. After the light beams have passed through the analysis zone, they are separated by a dichroic prism to be individually magnified and detected by one camera for each beam.

With more than 300 images per second and 4.2 megapixel per image the CAMSIZER X2 manages a 3.5 times higher data rate than the XT model while the software is capable of processing all images in real time. The new split view optical design is another highlight of the patented dual camera principle which is used in all CAMSIZER systems: One camera is optimized to

analyze the small particles with high resolution in a small field of view. The other camera detects the big particles in a large field of view, providing excellent statistics.

High precision reference object

Recalibrating the CAMSIZER X2 takes a matter of seconds with the help of a highly precise ($\pm 0,1 \mu\text{m}$) reticle, made by electron beam lithography which contains differently sized particles. Thus, the requirements of instrument calibration are fulfilled.



Particle Analyzer

CAMSIZER® X2


Modular Design for Optimized Measurement Conditions

The CAMSIZER X2's modular "X-Change" system offers three alternative dispersion methods, permitting the selection of the optimum method for each sample type:

- Air pressure dispersion
- Gravity dispersion
- Wet dispersion

The modules and plug-in cartridges are easily exchanged in one minute allowing for rapid switching between dispersion methods.

A video at www.retsche-technology.com/camsizer demonstrates how easily the "X-Change" system is operated.




Air Pressure Dispersion with "X-Jet"

The dispersion, that means the deagglomeration of the particles prior to passing through the measurement area, can be a crucial pre-condition for a proper measurement of the particles. Thanks to the flexible pressure adjustments of the "X-Jet" plug-in cartridge each material can be measured under optimum conditions.

Small particles have a relatively large surface area which increases the effects of adhesion forces or electrostatic charges, and can lead to aggregation of the particles. When passing through the nozzle, the sample material is dragged along by the compressed air and submitted to shear forces which break up agglomerates. These shear forces are enhanced by increasing pressure so that it is possible to separate even those particles which adhere tightly to each other. Too much pressure may be counter-productive: the shear forces may destroy the primary particles thus "grinding" the particles inside the nozzle.


With the dynamic image analysis method, it is possible to detect agglomerates or broken particles by analyzing the particle shape quantitatively or simply by looking at individual particle images. If comminution of the particles is happening, it will be necessary to adjust the pressure. Other measurement methods like laser diffraction use similar dispersion principles but are not capable of providing information on particle shape. The sample is collected in a vacuum cleaner after the measurement.

An optional cyclone is available for sample recovery. Although the air pressure dispersion accelerates the particles to up to 50 m/sec, the patented dual camera system ensures that wide particle size distributions as well as narrow, mono-modal samples below 10 microns are analyzed accurately.



Benefits at a glance

- Good for fine and agglomerated powders
- Pressure adjustable from 5 kPa to 450 kPa
- Automatic sample extraction with vacuum cleaner
- Optional sample recovery with cyclone



Particle Analyzer

CAMSIZER® X2

Gravity Dispersion with "X-Fall"

Flowable, unagglomerated samples can be analyzed by using the "X-Fall" plug-in cartridge. In this mode, the particles fall from a chute through the field of view of the two cameras accelerated only by gravity. Thanks to the comparably low speed of the particles, the large field of view and the high frame rate, the detection efficiency is very high, even for large particles up to 8 mm. Only a few coarser particles in the sample are sufficient for reliable, reproducible detection.

After the measurement, the sample material falls into a collection box and is available for additional analyses without loss or contamination.


Measuring range

10 μm  8 mm



Benefits at a glance

- Ideal for free flowing materials
- Complete sample recovery
- Gentle, contact-free measurement



Wet Dispersion with the "X-Flow" Module

The wet module "X-Flow" analyzes samples in the range from 0.8 µm to 1 mm in suspensions or emulsions. An advantage of this module is the small sample volume which is required. A low particle concentration in the dispersion medium, for example, 20 mg/l is already sufficient. It provides enough particles for a reproducible result in only 1 minute. Therefore, this module is highly suitable for applications involving pharmaceuticals or explosives. Here only a limited amount of sample material is available or advisable for safe operation. It can also be used for other "moist" samples including food and sand.

The measurement range of the "X-Flow" module starts at 0.8 µm. The CAMSIZER X2 also analyzes particles up to 1 mm without difficulty, provided they are kept suspended in the dispersion medium. Depending on the maximum particle size in

the sample, measurement cells of up to 4 mm may be used. Agglomerates can be further dispersed by an integrated ultrasonic probe.

For applications with organic solvents as the dispersion medium, the X-Flow module is equipped with a circulation system consisting of PTFE components and solvent-resistant seals.


Benefits at a glance

- Perfect for emulsions and suspensions
- Ideally suited for small sample volumes
- Resistant to organic solvents
- Integrated ultrasonic probe



Measuring range

0.8 µm — 1 mm






CAMSIZER® X2 Results

Evaluation and documentation

Retsch Technology supplies the CAMSIZER X2 with a powerful, process-oriented control and evaluation software which offers flexible data export to the customer's LIMS. A major advantage is the evaluation of the results in real-time. Graphical representation of the results is available while the measurement is still in progress. At the same time the measurement process can be checked visually by observing the digital images. All particle images and parameters are evaluated directly during measurement and are also saved by the software module "Particle X-Plorer" in a database. Immediately after the measurement is finished, the results can be displayed in various report forms. By using the particle database, it is possible to clarify particular issues in detail at a later date, and even to recalculate the results.

Simple and reliable operation

The user can easily select individual measurement and evaluation parameters and save product-specific settings. These are known as standard operating procedures (SOP), and they simplify the change between different repetitive measuring tasks. The SOP can be protected against corruption by a password protection. This ensures that the same instrument settings and output formats are always used with the highest degree of reliability and effectively eliminates operator-induced errors.




Examples of size parameters from the analysis of individual particle images derived from the high-resolution scans of the CAMSIZER X2

Particle shape analysis and its applications

In many application areas particle shape information is an important process and quality indicator. Based on Dynamic Image Analysis, the CAMSIZER X2 can simultaneously analyze the particle size and particle shape of the sample material in a sophisticated fashion.

Because of the high information content obtained from the digital images measured and stored, the particle projections can be evaluated in a variety of ways. Depending on the application, the CAMSIZER X2 measures various areas, perimeters and lengths of the particle projections, determining up to 50 different parameters per particle. Typical results obtained would be expressed as:

- Chord length
- Straight length
- Feret diameter
- Martin diameter
- Aspect ratio (width/length)
- Convexity
- Roundness
- Symmetry
- Transparency
- Angularity



In addition to different size definitions, various shape parameters are determined as significant indicators of the product quality. Those which are most important are illustrated in the graphics above.

Application Examples of Particle Shape Analysis:

- ▶ **Determination of edginess**
e. g. quality control of abrasives before further processing
- ▶ **Analysis of the broken fractions of granules**
e. g. replacement of time-consuming breakage analysis of the rolling properties of granules
- ▶ **Analysis of twins, triplets and agglomerates**
e. g. in pharmaceutical pellets
- ▶ **Determination of angularity**
e. g. geological analyses of sand or analyses of ceramic materials or abrasives
- ▶ **Analysis of fiber length and diameter**
e. g. of glass or carbon fibers
- ▶ **Determination of the roundness of metal powder particles**
e. g. quality control for 3D printing or of flowability and packing density of solder powder

A wealth of information from the advanced technology

To meet individual user requirements, the measurement results are displayed as graphs, tables, characteristics or particle images.

The results obtained can be presented graphically and in tabular form as size fractions, frequency distribution or cumulative distribution. In addition, the CAMSIZER X2 software allows the presentation of daily reports, trend analyses, mean value calculations and much more. A clearly structured, individually configurable measuring protocol based on international standards is produced for each analysis.


It is possible to store product specifications with regards to particle size and particle shape. If the measurement result deviates from these specifications, for example the roundness of the sample does not match the specifications, the software automatically produces a warning message.



Quality control

Comparison of the measurement result with upper and lower specification limits


Trend analysis of production processes
Up to 4 freely selectable sample parameters can be continuously monitored



NEW: Particle X-Plorer with 3D Cloud

The new structure of the CAMSIZER X2 software not only allows for data evaluation in real time but also for the **storage of large amounts of image data in the Particle X-Plorer database**. Particle images with noteworthy characteristics may be extracted and displayed at a later time. Similar image libraries are commonly used for documentation purposes in static image analysis (microscopy). This method allows for measurement of only a small number of particles. The CAMSIZER X2 particle database, however, handles millions of single particle images... for each measurement!

In addition to the well-known two-dimensional display of results, the **CAMSIZER X2 software supports three-dimensional graphics**, i. e. three different parameters being displayed on the same graphic (3D cloud). The cloud makes it possible to detect differences between samples which are not visible in a two-dimensional display. Moreover, particle groups with specific characteristics may be identified and separately evaluated. Thus, it is possible to recalculate, for example, the size distribution of all round and transparent particles in a sample.





3D Cloud: Graphic display of the 3 analysis parameters size, transparency, and aspect ratio for a mixture of road marking materials consisting of antiskid aggregates and glass beads

Electrodeposited copper



Water-atomized copper




All particle images and corresponding measurement parameters can be optionally stored in the Particle X-Plorer database. The entire measurement is additionally displayed in a diagram (right).

Key Features

1. Wide dynamic measuring range

Coffee powder: Thanks to the patented dual camera principle, the CAMSIZER X2 offers excellent resolution of fine particles. It also provides meaningful statistics which is particularly important for the detection of a small number of oversized particles. The graphic shows the comparison between three ground coffee samples which were ground quite differently. The red curve represents the narrow size distribution of a homogeneous sample while the green curve shows a much larger amount of over- and undersized grains. The blue curve is a mixture, showing a significant number of coarse particles larger than 1 mm.


An excessive amount of fines in the coffee powder could result in blocked filters. Too many coarse particles release less flavor into the beverage which makes the coffee taste less strong and aromatic.



2. High resolution of narrow distributions

Latex standards: The new patented split view optical design is particularly suitable for measuring narrow particle size distributions. The high resolution of the CAMSIZER X2 permits a very exact determination of the distribution breadths. Even in the lower measuring range of a few microns bimodal mixtures with very narrow size distributions are clearly resolved.



The example shows a mixture of two polystyrene latex standards with diameters of 2.5 μm and 5 μm respectively.



3. Outstanding repeatability and reproducibility

Metal powder: Even the smallest amounts of undersized or oversized particles are reliably detected, regardless of the measurement time. Moreover, the shape of small particles, for example the aspect ratio of particles smaller than 10 microns, is analyzed more precisely thanks to the high resolution of the CAMSIZER X2.


The example shows a fine metal powder with a size distribution from 5 μm to 20 μm like it is used, for example, in additive manufacturing, as solder powder or in other powder metallurgical manufacturing processes.



4. Particle shape analysis


Catalysts: Non-spherical (broken) particles can be detected by looking at parameters such as aspect ratio (breadth divided by length, b/l) or symmetry. The graphic shows a series of measurements with increasing amounts of broken particles during a feeding process, leading to an increase of the broken fraction in the aspect ratio diagram.

The initial material is almost spherical with a b/l ratio of more than 0.95 (red curve). The increasing fraction of broken particles can be derived directly from the graph at the threshold value 0.95. The actual particle images are shown on page 14, as example 9.



5. Size results 100% compatible with sieve analysis

Fluidized bed granules: Sieve analysis results may be compared directly with the data provided by CAMSIZER X2. The example shows the measurement of a granular mineral which is used as an additive to animal food. Thanks to the match between CAMSIZER X2 and sieve analysis results, product specifications can be directly compared even if they have been measured with different analytical techniques.




6. Comparison to laser diffraction and sieve analysis results

Granulated beverage powder: Dynamic image analysis simultaneously determines the particle size distribution with regards to width (red), length (green) and mean diameter (blue). The orange line represents the result of a measurement carried out with a laser diffraction analyzer. This corresponds roughly with the mean diameter obtained by dynamic image analysis.

Upon closer examination, the differences are apparent; the laser diffraction curve shows a broader distribution. For the large particles around 0.5 mm the laser diffraction measurement shows oversized particles which are not visible in the sieve analysis results. Measurement of the particle width is identical to sieve analysis. Laser diffraction measurements are based on the assumption that all particles are spherical. Therefore, it is not possible to differentiate parameters such as "diameter" and "length".


The curves shown in the example are typical for samples which contain non-round or longish particles.



Key Features

7. Reliable detection of oversized particles


PMMA micro beads: The CAMSIZER X2 features the latest camera technology with high frame rates and excellent resolution to capture as many particles as possible in the shortest time. The graphic shows the comparison of two samples with different amounts of oversized particles. Sample 2 (red) contains more oversized particles at 20 microns than sample 1 (green). The CAMSIZER's detection efficiency for small amounts of oversized particles is 500 times higher than that of any laser particle analyzer.



8. Dry measurement of small particles

Micro-grit abrasive: The X-Jet module is designed for the reproducible and exact measurement of fine powders. The example shows a sample of a SiC micro-grit abrasive at the lower end of the CAMSIZER X2 measuring range. The total size distribution lies between 1 micron and 12 microns. Thanks to a high depth of field, pulsed, ultra-strong light sources and short exposure times, even the rapidly flowing micron-sized particles are accurately detected.


The exactly defined fineness of the grinding material ensures an optimum combination of abrasion efficiency and surface roughness.



9. Analysis of individual particles: "What you see is what you get"

Catalyst spheres: With the Camsizer X2 software it is also possible to store and process images of individual particles. With one mouse click on the particle of interest all properties of this particle with respect to size and shape are displayed. Thus, the user can develop an intuitive understanding of the results, e. g. by learning which figures correspond to "good" or "bad" particles. This is especially important for method development to visually verify threshold values determined from the curves.



As an example, the software can detect agglomerated particles and uses them for the calculations or can also exclude them from the results.



10. Measurement of fiber length and diameter

Carbon fibers: The new optics with higher depth of field and resolution ensure consistently sharp images even for very thin, long fibers. Thus, it is possible to separately measure the fiber length and thickness quickly and reproducibly.

The example shows fibers with a diameter of less than 10 μm and a length of up to 2.5 mm. Fiber length and diameter are determined simultaneously thanks to the split view optical design.



Sieve Analysis, Laser Diffraction, Static or Dynamic Image Analysis?

A comparison of measurement techniques

Performance Features	CAMSIZER X2 Dynamic Image Analysis	Sieve Analysis	Laser Diffraction	opt. Microscope Static Image Analysis
1 Wide dynamic measurement range	++	+	++	-
2 Reproducibility and repeatability	++	+	++	-
3 High resolution for narrow distributions	++	-	-	++
4 Particle shape analysis	++	-	-	++
5 Direct measurement technique	++	++	-	++
6 Compatibility of results with other techniques	+	-	-	-
7 Reliable detection of oversized grains	+	++	-	-
8 Robust hardware, easy operation for routine analysis	++	++	++	-
9 Analysis of individual particles	+	-	-	++
10 High measurement speed, short measurement times	++	-	++	-

Specifications

Technical Data CAMSIZER® X2

Measurement ranges	Module "X-Dry" with "X-Fall" plug-in cartridge	10 µm to 8 mm
	Module "X-Dry" with "X-Jet" plug-in cartridge	0.8 µm to 5 mm
	Module „X-Flow“	0.8 µm to 1 mm
Measurement principle	Dynamic Digital Image Analysis (ISO 13322-2)	
Measurement time	approx. 1 to 3 min. (depends on required measurement statistics)	
Number of cameras	2 (measuring simultaneously)	
Sample volume	<20 mg – 500 g (depends on sample type and measurement mode)	
Measurement speed	>300 images/sec.	
Size of analysis area	approx. 20 x 20 mm	
Resolution	0.8 µm per Pixel	
Measurement parameters	Particle size	Smallest diameter, length, mean diameter etc.
	Particle shape	Aspect ratio (breadth to length), symmetry, sphericity, convexity etc. (ISO 9276-6)
Instrument data	Dimensions (H x W x D)	approx. 580 x 850 x 570 mm
	Weight (without PC)	approx. 50 kg
	Compressed air supply	approx. 6 - 8 bar
	Compressed air consumption	approx. 25 - 140 l/min (depending on dispersion pressure selected)
The CAMSIZER X2 is CE-certified and complies with the relevant guidelines and standards. The vacuum cleaner is included in the delivery scope of the X-Dry module.		
Options	Cyclone 21 CFR Part 11-compliant software IQ/OQ-documentation according to GLP/GMP	
Accessories	Different hoppers and chute materials allow for adaptation of the CAMSIZER X2 to the requirements of the particular sample material, depending on flowability and sample volume	




CAMSIZER® X2 at a glance

The CAMSIZER X2 analyzes the particle size and shape of fine powders, emulsions and suspensions with particle sizes of 0.8 µm and higher with excellent repeatability.

Easy operation, short measurement times and flexible dispersion options for agglomerated particles are standard. This allows for routine use not only in R&D but also in quality control laboratories with high sample throughput.

Compared to other particle size analysis methods, such as laser diffraction, the CAMSIZER X2 is characterized by its direct image processing measurement principle. This allows for analysis of representative sample volumes, even with wide size distributions in a very short time. The direct principle ensures a better understanding of the sample quality thanks to additional information (length, width, mean particle size and particle shape) and also a 500 times greater detection efficiency, e. g. for small quantities of oversized particles.



Retsch Technology GmbH
Retsch-Allee 1-5
42781 Haan, Germany

Phone +49 21 04 / 23 33 - 300
Fax +49 21 04 / 23 33 - 399

E-Mail technology@retschi.com
Internet www.retschi-technology.com

part of **VERDER**
scientific

RETSCH TECHNOLOGY – your specialist in particle analysis offers you a comprehensive range of instruments. We would be pleased to provide you with further information about our analytical instruments for dry and wet measurement with Dynamic Image Analysis or the Laser Diffraction method.