Quantachrome

Aquadyne DVS automated, gravimetric, water sorption analyzers

Quantachrome

- pharmaceuticals
- ► food
- construction materials
- fuel cell research
- papers
- carbons
- archaeological ceramics
- advanced materials



DVS-1/DVS-2





Aquadyne DVS-1/DVS-2 instrument features



1

Water sorption behavior plays an important role in the development and use of many common and advanced materials. Some examples include:

- Stability of pharmaceutical excipients and drug actives
- Drying and storage of grains
- Texture and shelf life of food products
- PEM fuel cell components
- Mortar, grout and other cementitious materials
- Paper and coatings
- Hydrophobic surface treatments
- Microporous and nanostructured carbons
- Archaeological ceramics

The Aquadyne DVS instruments are fully automated, gravimetric, one or two sample water vapor sorption analyzers. They measure adsorption and desorption isotherms of water vapor both accurately and sensitively, including sorption kinetics, with minimal operator involvement. The weight(s) of one or two sample(s) is / are constantly monitored and recorded as the relative humidity is automatically varied by the blending of dry carrier gas with a saturated gas stream using precision mass flow controllers. The dual balance design of the Aquadyne DVS-2 allows increased analysis throughput, side-by-side comparison with known or reference materials, or an extended mass range for a single sample.

The independently temperature controlled balance head environment ensures long term stability. The small sample chamber ensures rapid changes in sample atmosphere conditions when the relative humidity is altered during an analysis. Sample(s) can be pre-dried insitu up to 85° C (DVS-2) in a flow of dry gas. The software is available in a 21 CFR Part-11 compliant version, and runs under Windows[®] XP, Vista and Windows 7.

Measurement capabilities

Isotherms:

Mass change as a function of changing relative humidity, increasing mass during adsorption (increasing RH%), decreasing mass during subsequent desorption (decreasing RH%).

Kinetics:

Time-dependent studies give the rate of sorption.

Effect of Temperature:

Isotherms and kinetics change as a function of temperature. Can be used to yield sorption enthalpies.



Aquadyne DVS-1/DVS-2 applications/results

Deliquescence:

Certain salts exhibit the property of dissolving in self-adsorbed water.

Equilbrium moisture content:

Quantitatively the amount of water associated (adsorbed/ absorbed) with the material under given atmospheric conditions of temperature and relative humidity (organic materials may be subject to molding according to atmospheric humidity, for example).

Hysteresis working range:

Adsorption and desorption usually differ by hysteresis - this gives a range of stability in which subsequent sorption cycles produce no change in equilibrium moisture content.

Hydrophobicity/philicity:

The shape of the isotherm reveals relative strength of affinity between water and the surface, due to polar or other chemical interactions.

Micropore investigations:

Even hydrophobic materials can take up moisture if their pores are small enough.

Crystallization phenomena:

Certain amorphous materials will undergo crystallization as a function of water sorption due to, for example, shifts in glass transition points.

Sample form:

Sorption kinetics can be altered by the physical form of a material - powder, granule, pellet, tablet, monolith etc.



The dual sample capability of the Aquadyne DVS-2 is illustrated above.









Aquadyne DVS-1/DVS-2 software/specification

Software

The Aquadyne DVS is operated through userfriendly software running under Windows®. Sample analyses are easily and quickly set up in terms of target RH with equilibration criteria based on rate of weight change and/ or time. The software automatically records kinetic sorption data and isotherm points at user selectable intervals. During the run, data points yet to be recorded can be deleted, added to, or edited for different values (even sample chamber temperature) for real-time optimization.



Both kinetic and isotherm plots for points acquired can be viewed during the analysis, separately for each balance, or combined in an overlay plot. Once analyzed, results can be viewed in tabular and graphical format, and calculations such as BET surface area and heats of adsorption performed. Data can also be exported in a csv format.

Accessories

Regulator assembly, p/n: 01207

A part from electricity, the only utility the Aquadyne DVS requires is dry air or nitrogen regulated to the required pressure. Our two-stage regulator assembly includes CGA580 fitting and isolation valve.

Juantachrom



Specifications

DVS-1 Balance heads: 1 DVS-2 Balance heads: 2 Balances, type: ultrasensitive electronic microbalances Balance load capacity(each): 5g Balance load capacity(combined): 10g (DVS-2 only) Dynamic weighing range(each): -500mg to + 500mg, 0 to +1000mg Dynamic weighing range(combined): -1000mg to +1000mg, 0 to + 2000mg (DVS-2 only) Weighing resolution: 0.1µg RH range: <2% - 98% (depending on temperature). RH resolution: 0.1% RH accuracy: ± 0.8% RH at 25°C Sample chamber temperature range: <10 - 80°C (DVS-2), 15-60°C (DVS-1) Sample chamber temperature stability: ± 0.1°C Balance head temperature : Up to 85°C (DVS-2), 75°C (DVS-1) Gas flow rates: 0 - 200 cm³/min (mass flow controllers) Viewing port / sample chamber door: heated, triple-glazed Data acquisition rate: adjustable, up to 12 points per minute Sample pan material: silica (other materials can be employed) Water reservoir capacity: 46 ml

Utilities

Dimensions(W x D x H): 45cm x 60cm x 55cm (18" x 24" x 22") Weight: 60kg(132 lbs.)

Electrical: 100-230 V~, 50/60Hz - single phase Compressed gas: dry air or nitrogen regulated to 8-15 psig

PC: Windows 2000 or higher operating system with RS232 communications port and 9-pin D-sub connector.

Granite block, p/n: 90095 Aquadyne DVS base. Granite slab 46cm x 61cm x 8cm.



Cost 90 Project

In an effort to standardize the determination of water activity in foods, the European COST 90 project¹ made recommendations for the method of measurement of the water adsorption. The method specified using saturated salts to control the relative humidity and weighing the samples dry and after several days of equilibration.

Microcrystalline Cellulose (MCC) was selected as a standard and certified in a round-robin study involving ten independent labs². The mean percent adsorbed was determined and reported. These data form an internationally accepted standard for validation of instrumentation and methods used to measure water sorption.

Although these data were obtained using a static method, well-equilibrated dynamic vapor sorption (DVS) data should be in good agreement with these data.



302 data collected on an Aquadyne DVS-2 showing close agreement between the two balances. The black and red curves represent the CRM 302 mean adsorbed amount $\pm 2.77\sigma$.

Aquadyne DVS Performance

All Aquadyne DVS instruments are tested using the CRM 302 MCC standard and must fall within $\pm 2.77\sigma$ of the published mean before shipping. Validating these instruments to this internationally accepted standard ensures the highest quality data in your specific application.

CRM 302 data compared to acceptance limits

Data PointRH% ValueRH% LueSample TempMinActual AActual BMax111.1%11.1%25.001.7179632.0485381.9766552.542038222.5%22.5%25.002.7639763.2252223.1939483.716025333.0%33.0%25.003.8256334.109644.1147754.474367442.8%42.8%25.004.856274.8651234.9103165.463731552.9%52.9%25.005.4406535.6547975.7294256.499347657.7%57.7%25.005.9221226.1866016.2615787.037878770.8%75.3%25.107.9936188.7589118.9022869.806583984.3%84.4%25.109.79411210.7617210.9742511.838791090.2%90.1%25.0011.6001112.73912.9368414.93989								
111.1%11.1%25.001.7179632.0485381.9766552.542038222.5%22.5%25.002.7639763.2252223.1939483.716025333.0%33.0%25.003.8256334.109644.1147754.474367442.8%42.8%25.004.856274.8651234.9103165.463731552.9%52.9%25.005.4406535.6547975.7294256.499347657.7%57.7%25.005.9221226.1866016.2615787.037878770.8%70.8%25.007.6104078.0524478.175598.889593875.3%75.3%25.107.9936188.7589118.9022869.806583984.3%84.4%25.109.79411210.7617210.9742511.838791090.2%90.1%25.0011.6001112.73912.9368414.93989	Data Point	RH% Requested Value	RH% Actual- Value	Sample Temp	Min	Actual A	Actual B	Мах
222.5%22.5%25.002.7639763.2252223.1939483.716025333.0%33.0%25.003.8256334.109644.1147754.474367442.8%42.8%25.004.856274.8651234.9103165.463731552.9%52.9%25.005.4406535.6547975.7294256.499347657.7%57.7%25.005.9221226.1866016.2615787.037878770.8%70.8%25.007.6104078.0524478.175598.889593875.3%75.3%25.107.9936188.7589118.9022869.806583984.3%84.4%25.109.79411210.7617210.9742511.838791090.2%90.1%25.0011.6001112.73912.9368414.93989	1	11.1%	11.1%	25.00	1.717963	2.048538	1.976655	2.542038
3 33.0% 33.0% 25.00 3.825633 4.10964 4.114775 4.474367 4 42.8% 42.8% 25.00 4.85627 4.865123 4.910316 5.463731 5 52.9% 52.9% 25.00 5.440653 5.654797 5.729425 6.499347 6 57.7% 57.7% 25.00 5.922122 6.186601 6.261578 7.037878 7 70.8% 70.8% 25.00 7.610407 8.052447 8.17559 8.889593 8 75.3% 75.3% 25.10 7.993618 8.758911 8.902286 9.806583 9 84.3% 84.4% 25.10 9.794112 10.76172 10.97425 11.83879 10 90.2% 90.1% 25.00 11.60011 12.739 12.93684 14.93989	2	22.5%	22.5%	25.00	2.763976	3.225222	3.193948	3.716025
442.8%42.8%25.004.856274.8651234.9103165.463731552.9%52.9%25.005.4406535.6547975.7294256.499347657.7%57.7%25.005.9221226.1866016.2615787.037878770.8%70.8%25.007.6104078.0524478.175598.889593875.3%75.3%25.107.9936188.7589118.9022869.806583984.3%84.4%25.109.79411210.7617210.9742511.838791090.2%90.1%25.0011.6001112.73912.9368414.93989	3	33.0%	33.0%	25.00	3.825633	4.10964	4.114775	4.474367
552.9%52.9%25.005.4406535.6547975.7294256.499347657.7%57.7%25.005.9221226.1866016.2615787.037878770.8%70.8%25.007.6104078.0524478.175598.889593875.3%25.107.9936188.7589118.9022869.806583984.3%84.4%25.109.79411210.7617210.9742511.838791090.2%90.1%25.0011.6001112.73912.9368414.93989	4	42.8%	42.8%	25.00	4.85627	4.865123	4.910316	5.463731
657.7%57.7%25.005.9221226.1866016.2615787.037878770.8%70.8%25.007.6104078.0524478.175598.889593875.3%75.3%25.107.9936188.7589118.9022869.806583984.3%84.4%25.109.79411210.7617210.9742511.838791090.2%90.1%25.0011.6001112.73912.9368414.93989	5	52.9%	52.9%	25.00	5.440653	5.654797	5.729425	6.499347
770.8%70.8%25.007.6104078.0524478.175598.889593875.3%75.3%25.107.9936188.7589118.9022869.806583984.3%84.4%25.109.79411210.7617210.9742511.838791090.2%90.1%25.0011.6001112.73912.9368414.93989	6	57.7%	57.7%	25.00	5.922122	6.186601	6.261578	7.037878
8 75.3% 75.3% 25.10 7.993618 8.758911 8.902286 9.806583 9 84.3% 84.4% 25.10 9.794112 10.76172 10.97425 11.83879 10 90.2% 90.1% 25.00 11.60011 12.739 12.93684 14.93989	7	70.8%	70.8%	25.00	7.610407	8.052447	8.17559	8.889593
9 84.3% 84.4% 25.10 9.794112 10.76172 10.97425 11.83879 10 90.2% 90.1% 25.00 11.60011 12.739 12.93684 14.93989	8	75.3%	75.3%	25.10	7.993618	8.758911	8.902286	9.806583
10 90.2% 90.1% 25.00 11.60011 12.739 12.93684 14.93989	9	84.3%	84.4%	25.10	9.794112	10.76172	10.97425	11.83879
	10	90.2%	90.1%	25.00	11.60011	12.739	12.93684	14.93989

¹Spiess W.E.L. & Wolf W.R. (**1983**)" *The Results of the COST 90 Project on Water Activity in Physical Properties of Foods (Part I).*" Eds Jowett R., Eschen F., Halström B., Meffort H.F.Th., SPIES W.E.L. & VOS G. Applied Sci. Publishers, London.

²Jowitt, R. & Wagstaff, P.J. (**1989**)" *The Certification of the water content of microcrystalline cellulose (MCC) at 10 water activities — CRM 302.*" Office for Official Publications of the European Communities, Luxembourg.

4

CORPORATE HEADQUARTERS

Quantachrome Instruments

1900 Corporate Drive Boynton Beach, FL 33426 USA Phone: +1 (561) 731-4999 Fax: +1 (561) 732-9888 E-mail: qc.sales@quantachrome.com

www.quantachrome.com

EUROPE

Quantachrome UK Limited

 Pale Lane Farm, Pale Lane

 Hartley Wintney, Hook

 RG27 8BA, UK

 Phone:
 +44 (0)1252 819719

 Fax:
 +44 (0)1252 819901

www.quantachrome.co.uk





Quantachrome's Material Characterization Laboratory

EUROPE

Quantachrome GmbH & Co. KG

Rudolf-Diesel Str. 12 85235 Odelzhausen, Germany Phone: +49 (0)8134 93240 Fax: +49 (0)8134 932425 www.quantachrome.de

CHINA Quantachrome

Representative Office M806. Jingbao Garden

183 Andingmenwai Street, Beijing 100011, PRC Phone: +86 800 8100515

+86 10 64400892

+86 13 801191442

Fax: +86 10 64400892 www.quantachrome-china.com

JAPAN

· Gas Sorption Isotherms

True Solid Density

Tapped Density

Pore Size, Pore Size Distribution

Quantachrome is also recognized as

an excellent resource for authoritative

analysis of your samples in our fully equipped, state-of-the-art powder

(www.labqmc.quantachrome.com,

characterization laboratory, LabQMC

qc.lab@quantachrome.com).

Porosity, Pore Zeta Potential
Chemisorption, TPR/TPO/TPD
Water Sorption Behavior

Surface Area

Quantachrome Instruments Japan G.K.

KSP W311 Sakado 3-2-1, Takatsu-Ku, Kawasaki-shi, Kanagawa 213-0012, Japan Phone: +81 44 829 1525 Fax: +81 44 829 1527 www.quantachrome.co.jp

Quantachrome[®]

Renowned innovator of ideas for today's porous materials community.

For 45 years, Quantachrome's scientists and engineers have revolutionized measurement techniques and designed instrumentation to enable the accurate, precise, and reliable characterization of powdered and porous materials:



Quantachrome Instruments' quality management system is certified to be in accordance with ISO9001:2008.

Not only are Quantachrome products the instruments of choice in academia, but the technology conceived and developed by our expert staff is applied in industrial laboratories worldwide, where research and engineering of new and improved porous materials is ongoing. Manufacturers also rely on porous materials charac-terization technology to more precisely specify bulk materials, to control quality, and to isolate the source of pro-duction problems with greater efficiency.

Serving Porous Materials and Powder Characterization Needs Since 1968









www.quantachrome.com

Trademarks and registered trademarks are the property of their respective owners.

© 2013 Quantachrome Corporation 0313 07167 Rev A

-

Our policy of continuous development may cause the information and specifications contained herein to change without notice or liability.