

partica mini



Quality Control for Food • Pharmaceuticals • Cosmetics



Compact, Powerful Particle Size Analyzer

HORIBA's LA-350 Laser Diffraction Particle Size Distribution Analyzer is the ideal combination of performance, price, and packaging. Based on the advanced optical design of previous LA-series analyzers, the LA-350 strikes a harmonious balance between high-functionality, easy operation, low maintenance, and cost-effectiveness. The optimized design allows for a compact optical bench, resulting in an efficient use of bench space, while preserving the accuracy, precision and resolution that HORIBA's analyzers are famous for.



Powerful, convenient, and outstanding performance to meet your needs.




Small and powerful

The combination optical bench and circulation pump in one system is one of HORIBA's most popular designs. Now this design has a much smaller footprint which allows you to move the analyzer from place to place. This is especially valuable for quality control situations when the locations of sampling and analysis need to be separate to avoid contamination. Also, because it requires less space, it is possible to place the instrument where it is needed.


A powerful & versatile circulation system

The optical bench and circulation pump are combined into a single compact system. The compact size and low weight make this a convenient analyzer for today's crowded laboratories. It also makes it possible to transport easily to different locations within a facility or ship to remote locations for on-site testing where it is not practical to dedicate an analyzer.



A stable & reliable optical system

HORIBA's optical design ensures accuracy and stable measurement even in this small footprint. It features automatic alignment along with each blank. This prevents human error in the set up which allows even a novice user to obtain accurate measurements. The laser diode light source provides stable performance throughout the long lifetime of the analyzer. The dectectors, lens, and mirrors, which are the most sensitive components, are protected by placing them in the interior of the instrument. The design has been rigorously tested for durability and robustness.





Convenient to use

The software is simple to use. A single click selects the conditions and a second one generates the measurements. The system guides the user through the measurement, enabling stress free analysis without any advanced knowledge of the software. Additionally, there are flexible support features included in the standard software for those users who want to change the conditions and compare results with previous measurements.

$1\rightarrow 2\rightarrow 3$ step operation



Stress-free maintenance from an intelligent design

 The flow cell can be removed and replaced without any tools.



The flow cell cleaning is one of the most essential maintenance routine. Therefore it should be easily handled by anyone who operates the analyzer. The system is designed to keep enough space in the cell chamber for easy operations.

flow cell unit



The whole flow cell unit can be removed and replaced easily.

Small volume sampling or collecting samples?

Fraction cell unit



The Fraction cell optional accessory enables measurement of very limited sample amounts and collecting them after the measurement. The stir bar inside of the cell of the cell prevents the particle segregation and sedimentation.

• Optical alignment? Instant automatic alignment function with blanking and sample measurement ensures reproducible measurements and reliable performances.



Outstanding performance

Simple operation with exquisite performance!

The Partica mini covers a wide range of sizes: 100nm to 1000 microns.


The analysis guarantees that your production quality and development process will be accurate.

- Dynamic wide measurement range: 0.1-1000 microns
- Measurement accuracy support: ±1.4% guaranteed data accuracy with specified NIST traceable standard materials.
- ISO13320 compliant supporting criteria incorporated by HORIBA factory inspection procedure for each system.

| Specifications | |
|-----------------------------------|---|
| Measurement principle | Laser Diffraction and Mie Light Scattering Theory |
| Analysis materials | Powders, slurries, emulsions etc. |
| Measurement output | Particle diameters and size distribution, size related theoretical calculated values |
| Size measurement range | 0.1·1000μm |
| Analysis time | Typical measurement takes about 10 seconds from "Measure" to display the result. |
| Measurement method | Wet method: liquid dispersed particles with the flow sampling system |
| Required sample amount | 10mg-5g (depending on the sample size, distribution and materials) |
| Wet flow system liquid volume | approximately 130-230mL |
| Organic solvent compatibility | Available in solvent resistant flow sampling version or fraction cell system (Optional) |
| Measurement performance guarantee | HORIBA selected standard material use conditional support. |
| Operation | USB data communication with PC. |
| Data Processing/Results Display | Desktop or laptop PC/LCD; Printer |
| Operation conditions | 15°C∼35°C, 85% RH or less (Non-condensing) |
| Power | AC100/120/230V 50/60Hz, 150VA |
| Dimensions | W297mm×D420mm×H376mm (excluding the computer) |
| Mass | approximately 23kg |
| Optics | Light source: Laser diode 5mW, λ=650nm Analyzer Classification: Class 1A laser product Detectors: 64 ring detectors×1 Silicon photo detectors×6 |
| Wet sampling system | Ultrasonic: Ultrasonic probe inside of the flow system, 7 step power adjustment Circulation pumping system: Centrifugal pump, 15 steps of speed control Drain: Solenoid valve Flow cell material: Bolosilicate glass |
| Optional accessories | Fraction cell, Auto fill pump, Solvent resistant circulation system |

376 0 297 420

External Dimensions (mm)





Please read the operation manual before using this product to assure safe and proper handling of the product.

- The contents of this catalog are subject to change without prior notice, and without any subsequent liability to this company.
- The color of the actual products may differ from the color pictured in this catalog due to printing limitations. • It is strictly forbidden to copy the content of this catalog in part or in full.
- All brand names, product names and service names in this catalog are trademarks or registered trademarks of their respective companies.



http://www.horiba.com e-mail: info@horiba.co.jp

● HORIBA, Ltd. Head Office 2 Miyanohigashi, Kisshoin Minami-ku, Kyoto, Japan Phone: 81 (75) 313-8123 Fax: 81 (75) 321-5725

Tokyo Sales Office Kanda-Awaji-cho Nichome Building 2-6, Awaji-cho, Kanda, Chiyoda-ku, Tokyo, Japan

Phone: 81 (3) 6206-4721 Fax: 81 (3) 6206-4730

●HORIBA Instruments (Singapore) Pte Ltd. Head Office Hanoi Office 3 Changi Business Park Vista #01-01 Duy Tan Street, Duy Tan Street,

nanio Uffice Unit 10, 4 Floor, CMC tower, Duy Tan Street, Dich Vong Hau Ward, Cau Giay district, Hanoi, Vietnam Phone: 84 (4) 3795-8552 Fax: 84 (4) 3795-8553 Akzonobel House, Singapore 486051 Phone: 65 (6) 745-83 Fax: 65 (6) 745-8155

●HORIBA UK Ltd. 2 Dalston Gardens, Stanm Middx HA7 1BQ GREAT Stanmore. BRITAIN UK TEL:+44 (208) 204-8142 FAX: +44 (208) 204-6142

Bulletin:HRE-3684A

● HORIBA Jobin Yvon SAS 16-18 rue du Canal 91165, Longjumeau Cedex,

France Phone: 33 (1) 64-54-13-00 Fax: 33 (1) 69-09-07-21

●HORIBA (China) Trading Co., Ltd.
Shanghai Office
Unit D, 1F, Building A, Synnex
International Park, 1068
West Tianshan Road,
Shanghai, 200335 China
Phone: 86 (21) 6289-6060
Fax: 86 (21) 6289-5553

■Beijing Branch
12F, Metropolis Tower, No.2,
Beijing, 100080, China
Phone: 86 (10) 8567-9966
Fax: 86 (21) 6289-5553

● PT HORIBA Indonesia JI. Jalur Sutera Blok 20A, No. 16-17, Kel. Kunciran, Kec. Pinang Tangerang-15144, Indonesia Phone: 62 (21) 3044 8525 Fax: 62 (21) 3044 8521

● HORIBA Jobin Yvon GmbH GmbH Haupt-Str. 1* D-82008 Unterhaching, Germany Phone: 49 (89) 46-23-17-0 Fax: 49 (89) 46-23-17-99
 HORIBA Instruments Incorporated Irvine Office

 9755 Research Drive, Irvine, CA 92618, U.S.A.
 3880 Pa N. U.0882

 Phone: 1 (949) 250-4811
 Phone: Fax: 1 (749) 250-0924
 Edison Office 3880 Park Ave., Edison, NJ 08820, U.S.A. Phone: 1(732)494-8660 Fax: 1(732)549-5125

●HORIBA ITALIA SRL Via Luca Gaurico 209 - 00143, Roma, Italy Phone: 39 (6) 51-59-22-1 Fax: 39 (6) 51-96-43-34

● HORIBA Korea Ltd. 10, Dogok-Ro, 6-Gil, Gangnam-Gu, Seoul, 135-860, Korea Phone: 82 (2) 753-7911 Fax: 82 (2) 756-4972

●HORIBA India Private Ltd. Bangalore Office
Kamadhenu, No.17/1 - 32,
Bannerghatta Road,
Audugodi, Bangalore,
560030, India
Phone: 91 (80) 22210071

HORIBA

Printed in Japan 1510SK53

Guangzhou Branch Room 1611/1612, Goldion Digital Network Centar, 138 Tiyu Road East, Guangzhou, 510620 China Phone: 86 (20) 3878-1883 Fax: 86 (20) 3878-1810