



# HT3200A

## Top class GC autosampler

Establishing new standards for liquid autosamplers.



**HT3200A** has been engineered for **maximum performance and reliability**. The new system takes advantage of today's latest technology to deliver even greater reliability, performance, as well as offering an extraordinary GLP experience.

The **HT3200A** has a **top class sample capacity**: no other offers so many samples with near-to-**zero bench space** requirements. Due to its 209 2mL sample vials, the **HT3200A** is setting a new standard in the market.

Samples are organised in two racks that can be easily removed for sample loading or preparation,

or to be stored elsewhere. An **integrated Bar Code Reader** can be available with every **HT3200A** unit, offering full and detailed sample tracking.

The **HT3200A** features the innovative **SyringeID**, a proprietary technology<sup>1</sup> based on RFID tags. The SyringeID is able to **identify syringes** in an univocal way; thereby preventing errors when mounting a syringe, preventing syringe volume mismatching and keeps track of the syringe consumption (preventive maintenance counters). The SyringeID system is able to provide you with a level of confidence never previously

achieved by identification systems based on syringe carriers.

The **HT3200A** can be easily installed on all the GC and GC/MS systems available in the market, because of flexibility and modular configuration. It can be configured to **serve up to two injectors** in the most of the supported GCs.

The self-aligning "plug and play" **HT3200A** mounts in seconds, without tools. It can be easily moved between GCs when workloads change, due to its easy positioning, repositioning and easy removal.

## HTA - Sampling for science

[www.hta-it.com](http://www.hta-it.com)



**Easy to use:** Just load the samples and run the analysis with no extra downtime.

A large, full-color **touch screen interface** provides easier system accessibility and usability. The touch screen eliminates drilldown, simplifying instrument control for both novices and experienced users.

The **HT3200A is the fastest liquid autosampler available** on the market: injection is performed in less than 100ms. Fast-injection technology ensures the best possible peak shape, while maximising the accuracy of your results.

The **greater solvent capacity** means longer unattended operations. The autosampler supports the use of six solvent vials of 10ml each, giving a total capacity of 60ml.

Furthermore, advantage can be taken of the **double wash step capability**: as pre- and post-wash solvent in addition to A, B... F solvents you can also choose for a combination of A+B, A+C.... F+E solvents for superior analytical performance (carry-over adverse!).

The **HT3200A** also handles the **most sophisticated sampling techniques**, including the internal standard technique (also known as sandwich injection), multi-phase, ambient headspace, priority injection, nano-litre injection and much more. **Parameters are easily programmable** to optimise both the most convenient sampling methods for both extremely volatile or viscous samples and the best injection technique.

The **HT3200A** can help to process samples more quickly and to access better data. It can also mount **different types of syringes**, with volumes from 0.5 to 100µl.

**HT3200A** offers **syringe illumination** to always keep the sample under control, in order to check against air bubbles in method validations.

The **rotating tower leaves the injector port free for manual injection or maintenance**. The sample racks are mounted away from the GC oven to prevent exposure to high temperatures, which could cause degradation or condensation in the sample vial.

The **HT3200A** can be also controlled by a PC using the **HTA Autosampler Manager**.

## HT3200A Technical specifications

### General features

Syringe volume: 0.5, 1, 5, 10, 25, 50 and 100µl  
Tray capacity: 2 removable racks; 209 vials, 2ml  
Maintenance: preventive counters available  
Electrical control: LAN and TTL; optional: RS232  
Syringe area illumination: yes (programmable)  
SyringeID: included

### Filling

Sample volume: as low as step of 0.1µl  
Air volume: as low as step of 0.1µl  
Filling speed: 1-100µl/sec  
Viscosity delay: 0-15s  
Bubble elimination: up to 15 pull up strokes

### Injection

Injection speed: 1-100µl/sec  
Injection depth: programmable  
Pre and post inj delay: 0-99s

### Washing

Type: pre-injection, sample, post-injection  
Solvent capacity: 6x10ml vials  
Mode: single or double wash

### Internal standard technique

IS volume: as low as step of 0.1µl  
Air gap volume: as low as step of 0.1µl  
Mode: 1 or 2 air gaps

### Physical features

Dimensions (WxHxD)<sup>2</sup>: 280x570x320mm  
Weight: 9.2kg  
Power supply: 100-240±10%Vac; 50-60Hz; 60VA

<sup>2</sup> tray in closed position

<sup>1</sup> Patent pending

**HTA,**  
the  
company

HTA is one of the leading Italian scientific instrument engineering and manufacturing companies. We are currently focused on applications and solutions for analytical, life sciences and clinical chemistry automation. Our specialisation is in robotic systems for sample management; among our most popular products are the GC and HPLC autosamplers and preparative workstations. In addition, HTA offers engineering and manufacturing consultancy services for its OEM customers. HTA's quality management system is certified UNI EN ISO 9001:2008.



**Distributed by:**

**HTA s.r.l.**

via del Mella, 77-79 - 25131 Brescia - ITALY

T: +39 030 3582920 - F: +39 030 3582930

**info@hta-it.com - http://www.hta-it.com**