

Vibratory Sieve Shaker AS 450 basic

General Information

The AS 450 basic, is a budget-priced alternative to the AS 450 control sieve shaker. The new sieve shaker covers a size range from 25 µm to 125 mm and accepts loads of up to 15 kg. Time and amplitude are digitally set, a memory function allows storage of one program. The AS 450 basic is suitable for dry sieving. It is the economic solution for customers who need to sieve larger quantities of dry material with reliable results.

Application Examples

cement clinker, chemicals, coal, coke, construction materials, fillers, minerals, ores, plastics, sand, soils, ...

Product Advantages

- suitable for dry sieving
- excellent separation efficiency even with short sieving times
- efficient electromagnetic drive
- 3-D throwing motion which ensures optimum use of the open sieve area and lets the sample move equally over the whole sieving surface
- up to 10 fractions in one sieving operation
- free digital adjustment of all process parameters (time, amplitude, interval)
- mobile operation panel for comfortable handling
- maintenance-free
- easy operation, ergonomic design

Features

Applications	separation, fractioning, particle size determination
Field of application	chemistry / plastics, construction materials, engineering / electronics, environment / recycling, geology / metallurgy, glass / ceramics
Feed material	powders, bulk materials
Measuring range*	25 µm - 125 mm
Sieving motion	throwing motion with angular momentum - 3D movement
Max. batch / feed capacity	15 kg
Max. number of fractions	10 / 7
Max. mass of sieve stack	50 kg
Amplitude	digital, 0 - > 2 mm
Time display	digital, 0 - 99.9 min
Interval operation	1 - 99 s
Storable SOPs	1
Suitable for dry sieving	yes



Vibratory Sieve Shaker AS 450 basic

Suitable for wet sieving	-
Serial interface	-
Suitable sieve diameters	400 mm / 450 mm
Max. height of sieve stack	830 mm
Clamping devices	"standard"
Electrical supply data	different voltages
Power connection	1-phase
W x H x D	680 x 280 x 680 mm
Net weight	~ 140 kg
Standards	CE

*depending on feed material and instrument configuration/settings

Vibratory Sieve Shaker AS 450 basic

Videolink

Function Principle

The AS 450 basic works with an electromagnetic drive. This drive produces a 3D throwing motion that moves the product to be sieved equally over the whole sieving surface. The advantage: high stress capacity, extremely smooth operation and short sieving times with high separation efficiency.

Order data

Vibratory Sieve Shaker AS 450

for test sieves 400/450 mm (please order clamping device (AS 450 control), test sieves and collecting pan separately)

30.028.1001	AS 450 basic, 230 V, 50 Hz, incl. clamping device "standard"
-------------	--

Sieve stacks AS 450

60.166.000998	Sieve stack consisting of 10 test sieves (ISO 3310-1), 400 mm Ø, 65 mm height (63 µm, 125 µm, 250 µm, 500 µm, 1 mm, 2 mm, 4 mm, 8 mm, 16 mm, 31.5 mm) and collecting pan
60.168.000999	Sieve stack consisting of 7 test sieves (ISO 3310-1), 450 mm Ø, 100 mm height (63 µm, 125 µm, 250 µm, 500 µm, 1 mm, 2 mm, 4 mm) and collecting pan
60.167.000998	Sieve stack consisting of 10 test sieves (ASTM E11), 400 mm Ø, 65 mm height, (230 mesh, 120 mesh, 60 mesh, 35 mesh, 18 mesh, 10 mesh, 5 mesh, 5/16", 5/8", 1.1/4") and collecting pan
60.169.000999	Sieve stack consisting of 7 test sieves (ASTM E11), 450 mm Ø, 100 mm height (230 mesh, 120 mesh, 60 mesh, 35 mesh, 18 mesh, 10 mesh, 5 mesh) and collecting pan

Test Sieves 400 mm Ø

60.010.000400	Collecting pan, stainless steel, 400 mm Ø, height 65 mm
60.220.000400	Intermediate pan, stainless steel, 400 mm Ø, height 65 mm
60.935.000400	Intermediate ring, stainless steel, 400 mm Ø, height 65 mm
69.540.0050	Sieve lid, stainless steel, 400 mm Ø
05.114.0048	O-ring for test sieves, 400 mm Ø

Test Sieves 450 mm Ø

60.010.000450	Collecting pan, stainless steel, 450 mm Ø, height 50
---------------	--

Vibratory Sieve Shaker AS 450 basic

60.220.000450	mm Intermediate pan, stainless steel, 450 mm Ø, height 100 mm
60.935.000450	Intermediate ring, stainless steel, 450 mm Ø, height 100 mm
69.545.0050	Sieve lid, stainless steel, 450 mm Ø
05.114.0012	O-ring for test sieves, 450 mm Ø

400 mm Ø x 65 mm height

ISO 3310/1

	stainless steel with wire gauze
60.166.000025	25 µm *
60.166.000032	32 µm *
60.166.000036	36 µm
60.166.000038	38 µm
60.166.000040	40 µm *
60.166.000045	45 µm
60.166.000050	50 µm *
60.166.000053	53 µm
60.166.000056	56 µm
60.166.000063	63 µm *
60.166.000071	71 µm
60.166.000075	75 µm
60.166.000080	80 µm *
60.166.000090	90 µm
60.166.000100	100 µm *
60.166.000106	106 µm
60.166.000112	112 µm
60.166.000125	125 µm *
60.166.000140	140 µm
60.166.000150	150 µm
60.166.000160	160 µm *
60.166.000180	180 µm
60.166.000200	200 µm *
60.166.000212	212 µm
60.166.000224	224 µm
60.166.000250	250 µm *
60.166.000280	280 µm
60.166.000300	300 µm
60.166.000315	315 µm *
60.166.000355	355 µm
60.166.000400	400 µm *
60.166.000425	425 µm

Vibratory Sieve Shaker AS 450 basic

60.166.000450	450 µm
60.166.000500	500 µm *
60.166.000560	560 µm
60.166.000600	600 µm
60.166.000630	630 µm *
60.166.000710	710 µm
60.166.000800	800 µm *
60.166.000850	850 µm
60.166.000900	900 µm
60.166.001000	1.00 mm *
60.166.001120	1.12 mm
60.166.001180	1.18 mm
60.166.001250	1.25 mm *
60.166.001400	1.40 mm
60.166.001600	1.60 mm *
60.166.001700	1.70 mm
60.166.001800	1.80 mm
60.166.002000	2.00 mm *
60.166.002240	2.24 mm
60.166.002360	2.36 mm
60.166.002500	2.50 mm *
60.166.002800	2.80 mm
60.166.003150	3.15 mm *
60.166.003350	3.35 mm
60.166.003550	3.55 mm
60.166.004000	4.00 mm *
60.166.004500	4.50 mm
60.166.004750	4.75 mm
60.166.005000	5.00 mm *
60.166.005600	5.60 mm
60.166.006300	6.30 mm *
60.166.006700	6.70 mm
60.166.007100	7.10 mm
60.166.008000	8.00 mm *
60.166.009000	9.00 mm
60.166.009500	9.50 mm
60.166.010000	10.00 mm *
60.166.011200	11.20 mm
60.166.012500	12.50 mm *
60.166.013200	13.20 mm
60.166.014000	14.00 mm
60.166.016000	16.00 mm *

Vibratory Sieve Shaker AS 450 basic

60.166.018000	18.00 mm
60.166.019000	19.00 mm
60.166.020000	20.00 mm *
60.166.022400	22.40 mm
60.166.025000	25.00 mm *
60.166.026500	26.50 mm
60.166.028000	28.00 mm
60.166.031500	31.50 mm *
60.166.035500	35.50 mm
60.166.037500	37.50 mm
60.166.040000	40.00 mm *
60.166.045000	45.00 mm
60.166.050000	50.00 mm *
60.166.053000	53.00 mm
60.166.056000	56.00 mm
60.166.063000	63.00 mm *
60.166.071000	71.00 mm
60.166.075000	75.00 mm
60.166.080000	80.00 mm *
60.166.090000	90.00 mm
60.166.100000	100.00 mm *
60.166.106000	106.00 mm
60.166.112000	112.00 mm

125.00 mm *

450 mm Ø x 100 mm height

ISO 3310/1

	stainless steel with wire gauze
60.168.000036	36 µm
60.168.000038	38 µm
60.168.000040	40 µm *
60.168.000045	45 µm
60.168.000050	50 µm *
60.168.000053	53 µm
60.168.000056	56 µm
60.168.000063	63 µm *
60.168.000071	71 µm
60.168.000075	75 µm
60.168.000080	80 µm *
60.168.000090	90 µm

Vibratory Sieve Shaker AS 450 basic

60.168.000100	100 µm *
60.168.000106	106 µm
60.168.000112	112 µm
60.168.000125	125 µm *
60.168.000140	140 µm
60.168.000150	150 µm
60.168.000160	160 µm *
60.168.000180	180 µm
60.168.000200	200 µm *
60.168.000212	212 µm
60.168.000224	224 µm
60.168.000250	250 µm *
60.168.000280	280 µm
60.168.000300	300 µm
60.168.000315	315 µm *
60.168.000355	355 µm
60.168.000400	400 µm *
60.168.000425	425 µm
60.168.000450	450 µm
60.168.000500	500 µm *
60.168.000560	560 µm
60.168.000600	600 µm
60.168.000630	630 µm *
60.168.000710	710 µm
60.168.000800	800 µm *
60.168.000850	850 µm
60.168.000900	900 µm
60.168.001000	1.00 mm *
60.168.001120	1.12 mm
60.168.001180	1.18 mm
60.168.001250	1.25 mm *
60.168.001400	1.40 mm
60.168.001600	1.60 mm *
60.168.001700	1.70 mm
60.168.001800	1.80 mm
60.168.002000	2.00 mm *
60.168.002240	2.24 mm
60.168.002360	2.36 mm
60.168.002500	2.50 mm *
60.168.002800	2.80 mm
60.168.003150	3.15 mm *
60.168.003350	3.35 mm

Vibratory Sieve Shaker AS 450 basic

60.168.003550	3.55 mm
60.168.004000	4.00 mm *
60.168.004500	4.50 mm
60.168.004750	4.75 mm
60.168.005000	5.00 mm *
60.168.005600	5.60 mm
60.168.006300	6.30 mm *
60.168.006700	6.70 mm
60.168.007100	7.10 mm
60.168.008000	8.00 mm *
60.168.009000	9.00 mm
60.168.009500	9.50 mm
60.168.010000	10.00 mm *
60.168.011200	11.20 mm
60.168.012500	12.50 mm *
60.168.013200	13.20 mm
60.168.014000	14.00 mm
60.168.016000	16.00 mm *
60.168.018000	18.00 mm
60.168.019000	19.00 mm
60.168.020000	20.00 mm *
60.168.022400	22.40 mm
60.168.025000	25.00 mm *
60.168.026500	26.50 mm
60.168.028000	28.00 mm
60.168.031500	31.50 mm *
60.168.035500	35.50 mm
60.168.037500	37.50 mm
60.168.040000	40.00 mm *
60.168.045000	45.00 mm
60.168.050000	50.00 mm *
60.168.053000	53.00 mm
60.168.056000	56.00 mm
60.168.063000	63.00 mm *
60.168.071000	71.00 mm
60.168.075000	75.00 mm
60.168.080000	80.00 mm *
60.168.090000	90.00 mm
60.168.100000	100.00 mm *
60.168.106000	106.00 mm
60.168.112000	112.00 mm

Vibratory Sieve Shaker AS 450 basic

125.00 mm *

400 mm Ø x 65 mm height

ASTM E 11

	stainless steel with wire gauze
60.167.000025	25 µm; 500 Mesh
60.167.000032	32 µm; 450 Mesh
60.167.000038	38 µm; 400 Mesh
60.167.000045	45 µm; 325 Mesh
60.167.000053	53 µm; 270 Mesh
60.167.000063	63 µm; 230 Mesh
60.167.000075	75 µm; 200 Mesh
60.167.000090	90 µm; 170 Mesh
60.167.000106	106 µm; 140 Mesh
60.167.000125	125 µm; 120 Mesh
60.167.000150	150 µm; 100 Mesh
60.167.000180	180 µm; 80 Mesh
60.167.000212	212 µm; 70 Mesh
60.167.000250	250 µm; 60 Mesh
60.167.000300	300 µm; 50 Mesh
60.167.000355	355 µm; 45 Mesh
60.167.000425	425 µm; 40 Mesh
60.167.000500	500 µm; 35 Mesh
60.167.000600	600 µm; 30 Mesh
60.167.000710	710 µm; 25 Mesh
60.167.000850	850 µm; 20 Mesh
60.167.001000	1.00 mm; 18 Mesh
60.167.001180	1.18 mm; 16 Mesh
60.167.001400	1.40 mm; 14 Mesh
60.167.001700	1.70 mm; 12 Mesh
60.167.002000	2.00 mm; 10 Mesh
60.167.002360	2.36 mm; 8 Mesh
60.167.002800	2.80 mm; 7 Mesh
60.167.003350	3.35 mm; 6 Mesh
60.167.004000	4.00 mm; 5 Mesh
60.167.004750	4.75 mm; 4 Mesh
60.167.005600	5.60 mm; 3.5 Mesh
60.167.006300	6.30 mm; 1/4 "
60.167.006700	6.70 mm; 0.265 "
60.167.008000	8.00 mm; 5/16 "
60.167.009500	9.50 mm; 3/8 "

Vibratory Sieve Shaker AS 450 basic

60.167.011200	11.20 mm; 7/16 "
60.167.012500	12.50 mm; 1/2 "
60.167.013200	13.20 mm; 0.530 "
60.167.016000	16.00 mm; 5/8 "
60.167.019000	19.00 mm; 3/4 "
60.167.022400	22.40 mm; 7/8 "
60.167.025000	25.00 mm; 1 "
60.167.026500	26.50 mm; 1.06 "
60.167.031500	31.50 mm; 1 1/4 "
60.167.037500	37.50 mm; 1 1/2 "
60.167.045000	45.00 mm; 1 3/4 "
60.167.050000	50.00 mm; 2 "
60.167.053000	53.00 mm; 2.12 "
60.167.063000	63.00 mm; 2 1/2 "
60.167.075000	75.00 mm; 3 "
60.167.090000	90.00 mm; 3 1/2 "
60.167.100000	100.00 mm; 4 "
60.167.106000	106.00 mm; 4.24 "

125.00 mm; 5 "

450 mm Ø x 100 mm height

ASTM E 11

	stainless steel with wire gauze
60.169.000038	38 µm; 400 Mesh
60.169.000045	45 µm; 325 Mesh
60.169.000053	53 µm; 270 Mesh
60.169.000063	63 µm; 230 Mesh
60.169.000075	75 µm; 200 Mesh
60.169.000090	90 µm; 170 Mesh
60.169.000106	106 µm; 140 Mesh
60.169.000125	125 µm; 120 Mesh
60.169.000150	150 µm; 100 Mesh
60.169.000180	180 µm; 80 Mesh
60.169.000212	212 µm; 70 Mesh
60.169.000250	250 µm; 60 Mesh
60.169.000300	300 µm; 50 Mesh
60.169.000355	355 µm; 45 Mesh
60.169.000425	425 µm; 40 Mesh
60.169.000500	500 µm; 35 Mesh
60.169.000600	600 µm; 30 Mesh

Vibratory Sieve Shaker AS 450 basic

60.169.000710	710 µm; 25 Mesh
60.169.000850	850 µm; 20 Mesh
60.169.001000	1.00 mm; 18 Mesh
60.169.001180	1.18 mm; 16 Mesh
60.169.001400	1.40 mm; 14 Mesh
60.169.001700	1.70 mm; 12 Mesh
60.169.002000	2.00 mm; 10 Mesh
60.169.002360	2.36 mm; 8 Mesh
60.169.002800	2.80 mm; 7 Mesh
60.169.003350	3.35 mm; 6 Mesh
60.169.004000	4.00 mm; 5 Mesh
60.169.004750	4.75 mm; 4 Mesh
60.169.005600	5.60 mm; 3.5 Mesh
60.169.006300	6.30 mm; 1/4 "
60.169.006700	6.70 mm; 0.265 "
60.169.008000	8.00 mm; 5/16 "
60.169.009500	9.50 mm; 3/8 "
60.169.011200	11.20 mm; 7/16 "
60.169.012500	12.50 mm; 1/2 "
60.169.013200	13.20 mm; 0.530 "
60.169.016000	16.00 mm; 5/8 "
60.169.019000	19.00 mm; 3/4 "
60.169.022400	22.40 mm; 7/8 "
60.169.025000	25.00 mm; 1 "
60.169.026500	26.50 mm; 1.06 "
60.169.031500	31.50 mm; 1 1/4 "
60.169.037500	37.50 mm; 1 1/2 "
60.169.045000	45.00 mm; 1 3/4 "
60.169.050000	50.00 mm; 2 "
60.169.053000	53.00 mm; 2.12 "
60.169.063000	63.00 mm; 2 1/2 "
60.169.075000	75.00 mm; 3 "
60.169.090000	90.00 mm; 3 1/2 "
60.169.100000	100.00 mm; 4 "
60.169.106000	106.00 mm; 4.24 "

125.00 mm; 5 "