

General Information

The Rotor Beater Mill SR 200 is suitable for coarse and fine size reduction, either in batches or continuously. It can process dry, soft, medium-hard, organic and inorganic substances.

The SR 200 is intended for universal use: from sample preparation in laboratories up to preparing sample batches in pilot plants or production facilities.

Application Examples

bentonite, chemical products, coal, coke, drugs, feed pellets, fertilizers, food, graphite, gypsum, herbs, pharmaceutical materials, plant materials, plastics, powder coating, resins, rice, salts, seeds, soils, spices, sugar, ...

Product Advantages

- · suitable for batchwise operation of larger quantities
- high final fineness
- exchangeable grinding and sieve inserts
- optional grinding inserts 180° for grinding of hard-brittle materials by additional impact
- · optional distance rotor to reduce frictional heat
- quick-action door lock and motor brake
- defined final fineness due to bottom sieves with aperture sizes from 0.08
 10 mm
- maximum peripheral rotor speed 24.5 m/s

Features

Applications Field of application

Feed material Size reduction principle Material feed size* Final fineness* Speed at 50 Hz (60 Hz) Rotor peripheral speed Types of rotors Material of grinding tools

Sieve sizes

size reduction, deagglomeration agriculture, chemistry / plastics, construction materials, environment / recycling, food, geology / metallurgy, medicine / pharmaceuticals soft, medium-hard impact, shearing < 15 mm < 80 µm 2,850 min⁻¹ (3,420 min⁻¹) 20.4 m/s (24.5 m/s) rotor beater cast iron: grinding chamber stainless steel: rotor / ring sieves / retaining frame - chromium plated: retaining frame trapezoid holes 0.08 / 0.12 / 0.20 / 0.25 / 0.50 / 0.75 / 1.00 / 1.25 / 1.50





© Retsch GmbH - www.retsch.com - info@retsch.com Subject to technical modifications and errors page 1/5 07.03.2013



	/ 2.00 mm round holes 3.00 / 4.00 / 10.00 mm
Collector capacity	5 I / optional: 30 I
Drive	1-phase motor / 3-phase motor
Electrical supply data	different voltages
Power connection	1-phase / 3-phase
Protection code	IP 54
Power consumption	1.1 kW
W x H x D closed	560 x 1150 x 700 mm (with base frame)
Net weight	~ 57 kg (with base frame)
Documentation	Operation & Application Video
Standards	CE
*depending on feed material and instrument configuration/settings	



Videolink



http://www.retsch.com/sr200

Function Principle

Size reduction and deagglomeration in the SR 200 are achieved by impact and shear effects. The feed material passes from the hopper into the center of the grinding chamber where it is crushed between the rotor, sieve and grinding inserts.

As soon as the material is smaller than the aperture size of the sieve it enters the collecting receptacle. The quick-acting door lock ensures easy access to the grinding chamber for quick cleaning. With their motor brake (braking time 0.5 s), safety switch, splashback and access barrier in the inlet and outlet areas, the SR 200 offers the highest degree of operating safety.

Order data

Rotor Beater Mill SR 200 / SR 300

with rotor, 240 mm filter hose and collecting receptacle (please order retaining frame, sieves and optional base frame SR 200 separately)

20.732.0001	SR 200, $3/N$ ~ 400 V, 50 Hz, rotor speed 2,850 min ⁻¹ , of cast iron, without base frame
20.732.0003	SR 200, 230 V, 50 Hz, rotor speed 2,850 min ⁻¹ , of cast iron, without base frame

other electrical versions available for the same price

Retaining frames and ring sieves for SR 200 / SR 300

Retaining frame 360°	
22.642.0002	Retaining frame for ring sieves 360°, chromium plated
22.642.0001	Retaining frame for ring sieves 360°, of stainless steel

Ring sieves 360° of stainless steel

02.407.0057

Trapezoid holes, 0.08 mm

© Retsch GmbH - www.retsch.com - info@retsch.com Subject to technical modifications and errors page 3/5 07.03.2013



02.407.0026 02.407.0028 02.407.0030 02.407.0031 02.407.0032 02.407.0038 02.407.0033 02.407.0035 02.407.0036 02.407.0036

Retaining frame 180°

02.143.0014

Sieves 180° of stainless st	steel
-----------------------------	-------

03.647.0081 03.647.0039 03.647.0040 03.647.0041 03.647.0042 03.647.0043 03.647.0044 03.647.0045 03.647.0046 03.647.0047 03.647.0048 03.647.0049 03.647.0053

Additional items SR 200 / SR 300

01.824.0028	Base frame for SR 200
22.609.0003	Wheels for base frame
22.717.0003	Distance rotor of stainless steel for SR 200
22.717.0004	Distance rotor of stainless steel for SR 300
22.717.0001	Rotor of stainless steel for SR 200
22.717.0002	Rotor of stainless steel for SR 300
Additional items SR 200 / SR 300 / SK 100	

22.187.0001

Ring-type filter for collecting receptacle 5 litres

© Retsch GmbH - www.retsch.com - info@retsch.com Subject to technical modifications and errors page 4/5 07.03.2013

Trapezoid holes, 0.12 mm Trapezoid holes, 0.20 mm Trapezoid holes, 0.25 mm Trapezoid holes, 0.50 mm Trapezoid holes, 0.75 mm Trapezoid holes, 1.00 mm Trapezoid holes, 1.25 mm Trapezoid holes, 1.50 mm Trapezoid holes, 2.00 mm Round holes, 3.00 mm Round holes, 4.00 mm

Retaining frame with grinding insert for sieves 180°, stainless steel

Trapezoid holes, 0.08 mm		
Trapezoid holes, 0.12 mm		
Trapezoid holes, 0.20 mm		
Trapezoid holes, 0.25 mm		
Trapezoid holes, 0.50 mm		
Trapezoid holes, 0.75 mm		
Trapezoid holes, 1.00 mm		
Trapezoid holes, 1.25 mm		
Trapezoid holes, 1.50 mm		
Trapezoid holes, 2.00 mm		
Round holes, 3.00 mm		
Round holes, 4.00 mm		
Round holes, 10.00 mm, for pre-grinding		



22.748.0001	Dust filter clamping rings for ring-type filter, incl. 5 dust filters
22.524.0002	Dust filters for ring-type filter, 25 pieces
01.011.0023	Collecting receptacle, 5 litres, stainless steel
22.187.0003	Filter hose for collecting receptacle 5 litres, length 240 mm, incl. comfort flange
02.186.0027	Spare filter hose for 22.187.0003
22.003.0011	Collecting receptacle, 30 litres, plastic, incl. filter hose and comfort flange
02.186.0033	Spare filter hose for 22.003.0011