

#### **General Information**

The Rotor Beater Mill SR 300 is suitable for coarse and fine size reduction, either in batches or continuously. It can process dry, soft, medium-hard, organic and inorganic substances.

Its high speed in comparison to the SR 200 allows a significantly higher throughput and fineness.

The SR 300 is intended for universal use: from sample preparation in laboratories up to preparing sample batches in pilot plants or production facilities. The grinding chamber, the feed hopper and the material inlet are completely made from high quality stainless steel.

## **Application Examples**

bentonite, chemical products, coal, coke, drugs, feed pellets, fertilizers, food, graphite, gypsum, herbs, pharmaceutical materials, plant materials, plastics, powder coating, resins, rice, salts, seeds, soils, spices, sugar, ...

## **Product Advantages**

- · suitable for batchwise operation of larger quantities
- central locking device
- wide range of accessories including various collection and feeding systems, rotors and sieves
- high final fineness
- exchangeable grinding and sieve inserts
- grinding chamber, feed hopper and material inlet and outlet of stainless steel
- optional grinding inserts 180° for grinding of hard-brittle materials by additional impact
- optional distance rotor to reduce frictional heat
- · removable hopper for easy cleaning
- motor brake
- defined final fineness due to bottom sieves with aperture sizes from 0.08
  10 mm

## **Features**

Applications Field of application

Feed material Material feed size\* Final fineness\* Speed at 50 Hz (60 Hz) Rotor peripheral speed size reduction, deagglomeration agriculture, chemistry / plastics, construction materials, environment / recycling, food, geology / metallurgy, medicine / pharmaceuticals soft, medium-hard < 15 mm < 50 µm 8,100 min<sup>-1</sup> (9,700 min<sup>-1</sup>) 58 m/s (69.5 m/s)

© Retsch GmbH - www.retsch.com - info@retsch.com Subject to technical modifications and errors







Types of rotors	rotor beater / standard / distance	
Material of grinding tools	stainless steel: grinding chamber / ring sieves / retaining frame - chromium plated: retaining frame	
Sieve sizes	trapezoid holes 0.08 / 0.12 / 0.20 / 0.25 / 0.50 / 0.75 / 1.00 / 1.25 / 1.50 / 2.00 mm round holes 3.00 / 4.00 / 10.00 mm	
Collector capacity	5 I / optional: 30 I	
Drive	3-phase motor	
Electrical supply data	different voltages	
Power connection	3-phase	
Protection code	IP 54	
Power consumption	2.2 kW	
W x H x D closed	560 x 1200 x 890 mm (with base frame)	
Net weight	~ 95 kg	
Documentation	Operation & Application Video	
Standards	CE	
*depending on feed material and instrument configuration/settings		



#### Videolink



http://www.retsch.com/sr300

## **Function Principle**

Size reduction and deagglomeration in the SR 300 are achieved by impact and shear effects. The feed material passes from the hopper into the center of the grinding chamber where it is crushed between the rotor, sieve and grinding inserts.

As soon as the material is smaller than the aperture size of the sieve it enters the collecting receptacle. The quick-acting door lock ensures easy access to the grinding chamber for quick cleaning. With their motor brake (braking time < 0.5 s), safety switch, splashback and access barrier in the inlet and outlet areas, Rotor Beater Mills offer the highest degree of operating safety.

## Order data

#### Rotor Beater Mill SR 200 / SR 300

with rotor, 240 mm filter hose and collecting receptacle (please order retaining frame, sieves and optional base frame SR 200 separately)

20.733.1002

SR 300,  $3/N \sim 400 \text{ V}$ , 50 Hz, rotor speed 8,100 min<sup>-1</sup>, of stainless steel, with easily removable stainless steel hopper, incl. base frame

#### other electrical versions available for the same price

#### Retaining frames and ring sieves for SR 200 / SR 300

Retaining frame 360° 22.642.0002

22.642.0001

Retaining frame for ring sieves 360°, chromium plated Retaining frame for ring sieves 360°, of stainless steel

# **Ring sieves 360° of stainless steel** 02.407.0057

02.407.0057 02.407.0026

Trapezoid holes, 0.08 mm Trapezoid holes, 0.12 mm

page 3/5 07.03.2013



02.407.0028 02.407.0030 02.407.0031 02.407.0032 02.407.0038 02.407.0033 02.407.0034 02.407.0035 02.407.0036 02.407.0040

#### Retaining frame 180°

02.143.0014

#### Sieves 180° of stainless steel

03.647.0081 03.647.0039 03.647.0040 03.647.0041 03.647.0042 03.647.0043 03.647.0044 03.647.0045 03.647.0046 03.647.0048 03.647.0048 03.647.0049 03.647.0053

22.748.0001

#### Additional items SR 200 / SR 300

01.824.0028	Base frame for SR 200
22.609.0003	Wheels for base frame
22.717.0003	Distance rotor of stainless steel for SR 200
22.717.0004	Distance rotor of stainless steel for SR 300
22.717.0001	Rotor of stainless steel for SR 200
22.717.0002	Rotor of stainless steel for SR 300
Additional items SR 200 / SR 300 / S	K 100
22.187.0001	Ring-type filter for collecting receptacle 5 litres

Ring-type filter for collecting receptacle 5 litres Dust filter clamping rings for ring-type filter, incl. 5

© Retsch GmbH - www.retsch.com - info@retsch.com Subject to technical modifications and errors page 4/5 07.03.2013

Trapezoid holes, 0.20 mm Trapezoid holes, 0.25 mm Trapezoid holes, 0.50 mm Trapezoid holes, 0.75 mm Trapezoid holes, 1.00 mm Trapezoid holes, 1.25 mm Trapezoid holes, 1.50 mm Trapezoid holes, 2.00 mm Round holes, 3.00 mm Round holes, 4.00 mm

Retaining frame with grinding insert for sieves 180°, stainless steel

Trapezoid holes, 0.08 mm Trapezoid holes, 0.12 mm Trapezoid holes, 0.20 mm Trapezoid holes, 0.25 mm Trapezoid holes, 0.50 mm Trapezoid holes, 0.75 mm Trapezoid holes, 1.00 mm Trapezoid holes, 1.25 mm Trapezoid holes, 1.25 mm Trapezoid holes, 1.50 mm Trapezoid holes, 2.00 mm Round holes, 3.00 mm Round holes, 4.00 mm



	dust filters
22.524.0002	Dust filters for ring-type filter, 25 pieces
01.011.0023	Collecting receptacle, 5 litres, stainless steel
22.187.0003	Filter hose for collecting receptacle 5 litres, length 240 mm, incl. comfort flange
02.186.0027	Spare filter hose for 22.187.0003
22.003.0011	Collecting receptacle, 30 litres, plastic, incl. filter hose and comfort flange
02.186.0033	Spare filter hose for 22.003.0011

© Retsch GmbH - www.retsch.com - info@retsch.com Subject to technical modifications and errors