

1500

ShaQer™

Mixing and
Blending



Grinding and
Pulverizing



ShaQe up your QuEChERS method!

The ShaQer™ is the new, low cost solution for QuEChERS Sample Prep!

Benefits over handshaking include:

- » **Higher throughput** — shakes up to 6 samples simultaneously vs 2 by hand.
- » **Increased recovery of pesticide residue** in fruits and vegetables. The vigorous mixing action of the ShaQer™ results in thorough extraction.
- » **No variation in shaking process** — Ensures equal mixing conditions for all samples, which cannot be guaranteed when samples are shaken manually.
- » **Versatile**—Accommodates a variety of tube sizes and shapes.



See why the ShaQer™ is the ideal tool for the QuEChERS method!

Scan this code with the QR Reader on your smart phone or visit www.youtube.com/user/sampleprep



www.spexsampleprep.com
Call us at 1-855-GET-SPEX

SPEX® SamplePrep®

A Simple, Safe & Effective solution for QuEChERS



Safety interlock prevents unit from operating when lid is open. Window allows analyst to view samples during operation.

Vertical shaking motion ensures thorough mixing and analyte extraction.

Clamp holds six 50 mL vials or six 15 mL vials. Alternate clamp for six 50 mL vials or six 2 mL vials.

Spill tray for easy cleanup.

Digital timer display with adjustable operating time.

Compact, powerful motor shakes samples at 1500 rpm.



SPEX CertiPrep QuEChERS Kits

SPEX CertiPrep provides a complete solution to all your QuEChERS needs from a comprehensive set of Internal Standards for the QuEChERS analytical techniques (GC, GCMS, LCMS) to a catalog of over 4000 analytes that can be purchased as single components or mixes. Our customers can design their own custom standards with customizable kits that contain the materials necessary to take your samples all the way through the analytical process from preparation to testing.

For more information call 1-800-LAB-SPEX or visit www.spexcertiprep.com.



SPEX CertiPrep 

Corporate Headquarters

SPEX SamplePrep
15 Liberty Street
Metuchen, NJ 08840
1-855-GET-SPEX

Tel: 732-623-0465
Fax: 732-906-2492
sampleprep@spex.com
www.spexsampleprep.com

European Division

SPEX CertiPrep Ltd
2 Dalston Gardens
Stanmore, Middlesex HA7 1BQ

Tel: +44(0)208-204-6656
Fax: +44(0)298-204-6654
sales@spexcertiprep.co.uk
www.spexcertiprep.com

Your Science is Our Passion.®

